

Laboratory roll mixer

BS – P 06



- Tabular arrangement control panel allows easy options one's from 7 programme, indication justified parameter and operating reports..
- Top display showing time which was set /justified or subtraction time at minutes, bottom display showing you return turns of driving roll/min., or programme which you set. With membrane keyboard is possible to set programme 1-7.
- Adjustable distance between rolls allows apply at mixing different cylindrical bin from average 35 mm.
- Bin is possible apply in general height 500 mm (there are not a part of suppliers – it is possible to order them).

Technical data:

BS – P 06

Time	1 to 999 min (for P-7 dwell time and rotation 1 až 999 min)
Rotation speed of the roll	64 rpm
Programs	7 programmable combinations of total time, dwell time and partial rotation (parameter settings remain in memory of device even after switching off)
Weight	10,8 kg
Power	75 W
Power supply	230 V / 50 Hz
Fuse	1A / T35
Dimensions (l x w x h)	775 x 370 x 310 mm