

## LABORATORY HAMMER MILL

# PSY MP 20/ MP 40



The PSY performs a particle size distribution compatible with various type of analysis and products, and in accordance with several international norms (e.g. ICC and ISO for Falling time, Near Infrared and Gluten analysis).

### FEATURES

Hammer mill

Cyclone effect Self cleaning

Sieve 0.8mm

Compatible with NIR analysis

Compatible with Falling Time & Gluten analysis

Capacity 300 grams in 30-50 seconds

### SPECIFICATION

	<b>PSY MP 20</b>	<b>PSY MP 40</b>
Power supply	mono-phase 230 V / 50 Hz	tri-phase 400 V / 50 Hz
Motor power		1.5 kW
Rotation motor		2 840 rpm
Rotation hammer		12 000 rpm
Sieve (standard)		0.8mm
Size of driving belt		10A0710C
Diameter grinding room		104 mm
Depth grinding room		40 mm
Capacity		300 grams in less than 1 minute
Size inlet for granule		30 x 22 mm
Size inlet for powdery materials		30 x 10 mm
Volume of vassel for scraping sample		0,7 l
Dimensions		555 x 400 x 690 mm
Weight		39 Kg

### Producer:

MEZOS, spol. s r. o., Třída Edvarda Beneše 573, 500 12 Hradec Králové, Czech Republic  
tel/fax : (+420) 491 424 788

email: [mezos@mezos.cz](mailto:mezos@mezos.cz)

[www.mezos.cz](http://www.mezos.cz)